


Air Cooled Water Distiller

- DRYING OVEN
- LABORATORY INCUBATOR
- MUFFLE FURNACE
- WATER DISTILLER**
- WATER BATH & OIL BATH
- HEATING MANTLE
- STIRRER
- MIXER & SHAKER
- HEATING PLATE & SAND BATH
- DIGESTION APPARATUS & METAL BATH
- NITROGEN SAMPLE CONCENTRATOR
- CLEAN BENCH
- STERILIZER
- VACUUM PUMP
- FILTRATION
- GRINDER
- COLONY COUNTER
- DRYER
- CENTRIFUGE
- ELECTRIC FURNACE
- LIFTING TABLE
- LABORATORY CONSUMABLE



Dimension(mm)



DRYING OVEN

LABORATORY INCUBATOR

MUFFLE FURNACE

WATER DISTILLER

WATER BATH & OIL BATH

HEATING MANTLE

STIRRER

MIXER & SHAKER

HEATING PLATE & SAND BATH

DIGESTION APPARATUS & METAL BATH

NITROGEN SAMPLE CONCENTRATOR

CLEAN BENCH

STERILIZER

VACUUM PUMP

FILTRATION

GRINDER

COLONY COUNTER

DRYER


CENTRIFUGE

ELECTRIC FURNACE

LIFTING TABLE

LABORATORY CONSUMABLE

It is suitable for medical, health care, scientific research, chemical industry, colleges and laboratories to produce distilled water.



Colleges and universities



Life sciences




Fundamental research



Medical research



Industrial production



Epidemic prevention

Characteristics:

- Air-cooled shunt structure
- Large screen LCD display working status
- Automatic shutdown without water
- Distilled water can be produced by repeated preparation
- Easy to operate and reliable

Major Technical Parameters:

Model	SQF-4
Output Quantity (L/H)	4
Power (KW)	4.5
Voltage	110V/220V
Overall Dimension (L*W*H)(mm)	490*490*1170
Weight (kg)	61 KG
Safety Device	Automatic power-off for water shortage
Additional Function	Air-cooled shunt structure
Type	Air-cooled
Cooling fan	Low noise axial flow (low noise outer rotor fan) Air