

# Desktop Laboratory Incubator

DRYING OVEN

LABORATORY INCUBATOR

MUFFLE FURNACE

WATER DISTILLER

WATER BATH & OIL BATH

HEATING MANTLE

STIRRER

MIXER & SHAKER

HEATING PLATE & SAND BATH

DIGESTION APPARATUS & METAL BATH

NITROGEN SAMPLE CONCENTRATOR

CLEAN BENCH

STERILIZER

VACUUM PUMP

FILTRATION

GRINDER

COLONY COUNTER

DRYER

CENTRIFUGE


ELECTRIC FURNACE

LIFTING TABLE

LABORATORY CONSUMABLE



## Dimension(mm)



DRYING OVEN

LABORATORY INCUBATOR

MUFFLE FURNACE

WATER DISTILLER

WATER BATH & OIL BATH

HEATING MANTLE

STIRRER

MIXER & SHAKER

HEATING PLATE & SAND BATH

DIGESTION APPARATUS & METAL BATH

NITROGEN SAMPLE CONCENTRATOR

CLEAN BENCH

STERILIZER

VACUUM PUMP

FILTRATION

GRINDER

COLONY COUNTER

DRYER

CENTRIFUGE

ELECTRIC FURNACE

LIFTING TABLE

LABORATORY CONSUMABLE

**Desktop Laboratory Incubator is suitable for the cultivation of the bacteria, microbial and scientific research, is the modern medicine, medicine, biological and biochemical industries test of the important equipment.**


Vaccine

Medicine

Environmental Science

Clinical Medicine

**Temperature Fluctuation Sheet:**



DRYING OVEN

### Characteristics:

LABORATORY INCUBATOR

- Novel vertical design, high quality cold-rolled steel plate production, electrostatic spraying process on the surface, modelling beautiful, strong and durable.
- With the latest European technology design, stainless steel tank, the corners of a circular arc transition, adjustable shelf, easy to clean.
- The working chamber is made of high quality stainless steel or cold rolling steel plate after anticorrosive processing, processing and into.
- With new structure, the door is in the middle is equipped with double tempered glass large Angle observation window, with elastic adjustable door lock, lock type high temperature sealing strip, ensure good sealing, prevent heat loss.
- In the workroom, the hot air is convective naturally and cyclically.
- Digital display microcomputer intelligent PID temperature controller.
- New anti-hot handle.

MUFFLE FURNACE

WATER DISTILLER

WATER BATH & OIL BATH

HEATING MANTLE

### Major Technical Parameters:

STIRRER

MIXER & SHAKER

HEATING PLATE & SAND BATH

DIGESTION APPARATUS & METAL BATH

NITROGEN SAMPLE CONCENTRATOR

CLEAN BENCH

STERILIZER

VACUUM PUMP

FILTRATION

GRINDER

COLONY COUNTER

DRYER

CENTRIFUGE

ELECTRIC FURNACE

LIFTING TABLE

LABORATORY CONSUMABLE

Model	WP-25A	WP-30A
	WP-25AB	WP-30AB
Cycle Mode	Natural convection	
Voltage	110/220V	
Volume (L)	18	30
Power(W)	200	300
Temp. Range (°C)	RT+5~ 65°C	
Temp. Accuracy (°C)	0.1°C	
Temp. Fluctuation (°C)	±0.5°C	
Temp. Uniformity (%)	≤4.7	
Timing Range	1-9999min/1 ~ 9999h (Constnt temperature time)	
Sensor	CU50	
Operation Function	Auto tuning function, Offset value setting, Timing function;	
Working Room Size (W*D*Hmm)	260*260*260	315*315*315
Overall Dimensions (W*D*Hmm)	370*420*550	420*480*650
Size of Package (W*D*H)(mm)	620*440*485	530*470*690
Number Of Shelves (pc)	2	2
Shelf Load (kg)	15	15
Max Shelf number (Standard 2pcs)	6	3
Spacing Of Shelf (mm)	35	80
NW/GW(kg)	16/18	22/25
Inner Chamber	A model: High strength rolled galvanized steel inner chamber; AB model: Mirror stainless steel inner chamber	
Outer Shell	Cold rolling steel electrostatic spraying exterior	
Insulation Layer	High quality rock wool board	
Heater	Ferrochrome wire	
Exhaust Hole diameter	36mm	
Safety Device	Over temperature sound-light alarm	
Optional Accessories	LCD controller, LCD programmable controller, Independent temperature limiter	