

Electric Furnace/Closed Electric Furnace

- DRYING OVEN
- LABORATORY INCUBATOR
- MUFFLE FURNACE
- WATER DISTILLER
- WATER BATH & OIL BATH
- HEATING MANTLE
- STIRRER
- MIXER & SHAKER
- HEATING PLATE & SAND BATH
- DIGESTION APPARATUS & METAL BATH
- NITROGEN SAMPLE CONCENTRATOR
- CLEAN BENCH
- STERILIZER
- VACUUM PUMP
- FILTRATION
- GRINDER
- COLONY COUNTER
- DRYER
- CENTRIFUGE
- ELECTRIC FURNACE**
- LIFTING TABLE
- LABORATORY CONSUMABLE



FL-1



FL-2



DL-1



It is suitable for heating in industry, agriculture, industry and mining, colleges and universities, families, medical and health care, scientific research units and laboratories.



Chemical Industry



University



Research Laboratories



Medical Industry



Laboratory

DRYING OVEN

LABORATORY INCUBATOR

MUFFLE FURNACE

WATER DISTILLER

WATER BATH & OIL BATH

HEATING MANTLE

STIRRER

MIXER & SHAKER

HEATING PLATE & SAND BATH

DIGESTION APPARATUS & METAL BATH

NITROGEN SAMPLE CONCENTRATOR

CLEAN BENCH

STERILIZER

VACUUM PUMP

FILTRATION

GRINDER

COLONY COUNTER

DRYER

CENTRIFUGE

ELECTRIC FURNACE

LIFTING TABLE

LABORATORY CONSUMABLE

Characteristics:

- The shell is made of high-quality cold-rolled steel plate, and the surface is treated by spraying process, which has strong corrosion resistance and is durable.
- Electronic temperature control device, with stepless temperature control, to meet the needs of different heating temperatures.
- The closed type is equipped with cast iron cover, which has the characteristics of heating without open flame.

Major Technical Parameters:

Name	Model	Specification	Power (w)	Overall dimension (l*w*h)(mm)	Nw (kg)	Gw (kg)
Electric furnace	DL-1	Single	1000	162*180*130	2.3	2.7
Electric furnace	DL-1	Single	2000	182*202*130	2.7	2.8
Electric furnace	DL-1	Two units	2000	300*180*130	4	4.4
Electric furnace	DL-1	Four units	4000	576*180*130	7.4	7.7
Electric furnace	DL-1	Six units	6000	855*180*130	10.9	11.3
Closed electric furnace	FL-1	Single	500	150*170*135	1.9	2
Closed electric furnace	FL-2	Single	1500	290*270*80	2.1	2.3