

Transparent water bath

DRYING OVEN

LABORATORY INCUBATOR

MUFFLE FURNACE

WATER DISTILLER

WATER BATH & OIL BATH

HEATING MANTLE

STIRRER

MIXER & SHAKER

HEATING PLATE & SAND BATH

DIGESTION APPARATUS & METAL BATH

NITROGEN SAMPLE CONCENTRATOR

CLEAN BENCH

STERILIZER

VACUUM PUMP

FILTRATION

GRINDER

COLONY COUNTER


DRYER

CENTRIFUGE


ELECTRIC FURNACE

LIFTING TABLE

LABORATORY CONSUMABLE



Dimension(mm)



DRYING OVEN

LABORATORY INCUBATOR

MUFFLE FURNACE

WATER DISTILLER

WATER BATH & OIL BATH

HEATING MANTLE

STIRRER

MIXER & SHAKER

HEATING PLATE & SAND BATH

DIGESTION APPARATUS & METAL BATH

NITROGEN SAMPLE CONCENTRATOR

CLEAN BENCH

STERILIZER

VACUUM PUMP

FILTRATION

GRINDER

COLONY COUNTER

DRYER

CENTRIFUGE

ELECTRIC FURNACE

LIFTING TABLE

LABORATORY CONSUMABLE

This equipment applies to medical unit, college, scientific research institution and the department of chemical printing, dyeing, pharmacy, petroleum and electronic component production for the use of evaporation, a refraction, concentration and constant temperature heating.




A must for college laboratories

Necessary for mining enterprises to experiment

Essential for medical research

Temperature Fluctuation Sheet:



Characteristics:

- Using innovative design structure, effectively prevent steam leakage.
- Intelligent temperature control instrument with high brightness digital display, light-touch type key has the advantages of simple operation and high precision temperature control.
- The whole structure is made of transparent engineering plastic, can clearly see the heated object.

Major Technical Parameters:

Model	TMS-2	TMS-4	TMS-5	TMS-7	TMS-9	TMS-12	TMS-20	TMS-26
Range of temperature control(°C)	RT+10~70							
Temp. Resolution Ratio	0.1°C							
Temp. Fluctuation (°C)	±1							
Temp. Accuracy	≤±1°C							
Voltage	110/220V							
Tming range	1min ~ 999 h							
Volume (L)	2	4	5	7	9	12	20	26
Power (W)	300	500	500	500	1000	1000	1500	1500
Size of Working Room (L*W*H)(mm)	175*160*97	265*160*97	325*175*97	325*175*150	325*265*97	325*265*150	530*325*97	530*325*150
Overall dimension (L*W*H)(mm)	275*230*165	358*230*165	420* 230*180	420*275*180	420* 230*270	420*275*270	625*230*330	627*275*330
Safety device	Effectively prevent steam leakage							
Optional accessories	LCD Display, Computer interface and printer.							